Monoclonal Mouse
Anti-Human
Cytokeratin 18
Clone DC 10
Code No./ Code/ Code-Nr. M 7010
Edition/ Ausgabe 27.02.03

ENGLISH

Intended use

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human Cytokeratin 18, Clone DC 10, is intended for use in immunocytochemistry. The antibody labels epithelial cells in those epithelia expressing cytokeratin 18, and hence the antibody is a useful tool for the differentiation and identification of epithelial tumours (1) and epithelioid hemangioendotheliomas (2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Introduction

Cytokeratin 18 belongs to the intermediate filaments, which create a cytoskeleton in almost all cells. In contrast to other intermediate filaments, cytokeratins (CKs) are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. Until now, 20 distinct CK polypeptides (3, 4) have been revealed in various human epithelial cells (5). They can be divided into an acidic (type I) and a neutral-basic (type II) subfamily. CK 18, a 45 kDa protein, belongs to the acidic type of cytokeratins, and is typically expressed in simple, nonstratified epithelia. However, CK 18 is also expressed in basal and superficial cells of transitional epithelium as well as in the luminal/secretory cells of complex epithelia (5).

Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN₃.

Clone: DC 10 (6). Isotype: IgG1, kappa.

Mouse IgG concentration: See label on vial.

Immunogen

PCM 42, human breast carcinoma cell line (6).

Specificity

In Western blotting of cytokeratin preparations from BT 20 cells, the antibody labels a single band of 45 kDa corresponding to CK 18 (6). In immunocytochemistry on human cultured cell lines, the antibody labels tumour cells of epithelial origin, whereas normal fibroblasts and tumour cells of non-epithelial origin such as glioma, melanoma and osteosarcoma, are not labelled (6).

As demonstrated by immunocytochemistry, the antibody cross-reacts with the CK 18-equivalent protein in monkey (7). No cross-reactivity with the CK 18-equivalent protein is observed in cow, dog, hamster, mouse, rat, sheep and swine (6).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Specimen preparation

Paraffin sections: The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or methacarn (6). Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval of tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, or 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labelling frozen sections (6).

Staining procedure

Dilution: Monoclonal Mouse Anti-Human Cytokeratin 18, code No. M 7010, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human kidney or mamma carcinoma and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen.

Visualization: DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

Automation: The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.

Performance characteristics

Cells labelled by the antibody display a cytoplasmic staining pattern.

Normal tissues: In general, the antibody labels simple epithelia and transitional epithelium as opposed to stratified epithelial components, which are negative. The labelling includes simple epithelia from the intestine, bronchial tree and alveoli of the respiratory system, kidney tubuli, bile ducts, surface mucosa of endometrium and endocervix, Fallopian tube epithelium, and rete epithelium of testis and ovary. Furthermore, hepatocytes and pancreatic acini are labelled by the antibody. In addition, the antibody labels the complex glandular epithelia of breast-, salivary-, and sweat glands, as well as the transitional epithelium of the urinary bladder (6). Also the epithelial cells of the periodontium are labelled (7). Of non-epithelial cells, the antibody labels endothelial cells in venules, lymphatics and capillaries in such

tissues as skin, subcutaneous soft tissues, skeletal muscle, placenta, and, e.g. respiratory, gastrointestinal, and genital tract mucosa (2). No labelling is observed in epithelial tissues which are known to lack CK 18, e.g. epidermis, foot sole epidermis, vagina, exocervix, esophagus, and myoepithelial cells (6).

Abnormal tissues: In tumours of pancreas, the antibody labelled 6/6 acinar cell carcinomas, 1/1 pancreatoblastoma, 19/20 pancreatic neuroendocrine tumours, and 10/19 solid-pseudopapillary tumours (8). Furthermore, the antibody labelled 30/30 biphasic, 21/46 monophasic- and 8/17 poorly differentiated synovial sarcomas (9). In vascular tumours, the antibody labelled endothelial cells in 17/17 epithelioid hemangioendotheliomas, 9/14 epithelioid angiosarcomas, 6/10 spindle-cell hemangiomas, 2/7 lymphangiomas, 10/48 non-epithelioid angiosarcomas, 3/13 venous hemangiomas, and 1/18 capillary hemangiomas. No labelling was observed in 4 cavernous hemangiomas and 6 Kaposi's sarcomas (2).

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1. Lane EB, Alexander CM. Use of keratin antibodies in tumor diagnosis [review]. Semin Cancer Biol 1990;1:165-79.
2. Miettinen M, Fetsch JF. Distribution of keratins in normal endothelial cells and a spectrum of vascular tumors: Implications in tumor diagnosis. Hum Pathol 2000;31:1062-7.
3. Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: Patterns of expression in normal epithelia, tumors and cultured cells [review]. Cell 1982;31:11-24.
4. Moll R, Löwe A, Laufer J, Franke WW. Cytokeratin 20 in human carcinomas. A new histodiagnostic marker detected by monoclonal antibodies. Am J Pathol 1992;140:427-47.
5. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors [review]. Subcell Biochem 1998;31:205-62.
6. Laurová L, Kovářík J, Bártek J, Rejthar A, Vojtěšek B. Novel monoclonal antibodies defining epitope of human cytokeratin 18 molecule. Hybridoma 1988;7:495-504.
7. Sculean A, Berakdar M, Pahl S, Windisch P, Brex M, Reich E, et al. Patterns of cytokeratin expression in monkey and human periodontium following regenerative and conventional periodontal surgery. J Periodont Res 2001;36:260-8.
8. Notohara K, Hamazaki S, Tsukayama C, Nakamoto S, Kawabata K, Mizobuchi K, et al. Solid-pseudopapillary tumor of the pancreas: immunohistochemical localization of neuroendocrine markers and CD10. Am J Surg Pathol 2000;24:1361-71.
9. Miettinen M, Limon J, Niezabitowski A, Lasota J. Patterns of keratin polypeptides in 110 biphasic, monophasic, and poorly differentiated synovial sarcomas. Virchows Arch 2000;437:275-83.

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