Monoclonal Mouse Anti-Human Cytokeratin 20
Clone K_s 20.8
Code M7019

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Cytokeratin 20, Clone K _s 20.8, is intended for use in immunohistochemistry. Cytokeratin 20 (CK 20) expression is almost entirely confined to normal and abnormal gastric and intestinal epithelium, urothelium and Merkel cells (1, 2), and the antibody is of particular value in the classification of carcinomas originating from these cell types (1, 3, 4, 5). Differential identification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	Protein IT (6).
Summary and explanation	CK 20 belongs to the intermediate filament proteins, which create a cytoskeleton in almost all cells. In contrast to other intermediate filaments, cytokeratins (CKs) are made up of a highly complex multigene family of 40 to 68 kDa polypeptides (1, 7). 20 distinct CK polypeptides have been revealed in various human epithelial cells and their malignant counterparts (1, 2, 7). They can be divided into an acidic (type I) and a neutral-basic (type II) subfamily. CK 20, a 46 kDa protein, is less acidic than the other type I cytokeratins. CK 20 positivity (5-100% CK-20-positive tumour cells) has been reported in the vast majority of adenocarcinomas of the colon, mucinous ovarian tumours, transitional-cell and Merkel cell carcinomas, and frequently also in adenocarcinomas of the stomach (50%), bile system (75%), and pancreas (40%). Less than 2% of breast and lung adenocarcinomas and squamous cell carcinomas were positive, while adenocarcinomas of the ovary (non-mucinous) and endometrium, and renal cell and small cell carcinomas showed only scattered positive cells in some cases (1).
Reagent provided	Monoclonal mouse antibody provided in liquid form as purified mouse IgG from ascitic fluid. In 0.05 mol/L Tris/HCL, 15 mmol/L NaN ₃ , 1 % bovine serum albumin, pH 7.2. <u>Clone:</u> K _s 20.8 (1). <u>Isotype:</u> IgG2a, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
Immunogen	CK 20 isolated by SDS-PAGE from cytoskeletal material obtained from human duodenal mucosa (1).
Specificity	Anti-Human CK 20, Clone K _s 20.8, labels a 46 kDa polypeptide corresponding to CK 20 in Western immunoblotting of SDS-PAGE-separated cytoskeletal proteins from human duodenal mucosa (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling formalin-fixed, paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with trypsin (4), proteinase K or heat-induced epitope retrieval is required. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, High pH, code S3308, Dako Target Retrieval Solution, code S1700, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody cannot be used for labelling frozen sections or cell smears owing to cross-reactivity with some cytokeratin-20-negative epithelia (e.g. renal tubules) (1).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Cytokeratin 20, code M7019, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of human colon, and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, High pH, code S3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG2a, code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> Dako LSAB™+/HRP kit, code K0679, and Dako EnVision™+/HRP kits, codes K4004 and K4006, are recommended. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunohistochemical staining using automated platforms, such as the Dako Autostainer.
Product-specific limitations	CK 20 may occasionally be expressed in breast and lung adenocarcinomas, and in squamous cell carcinomas (1, 3). Less than 5% CK-20 positive cells may be present in a number of tissues not generally considered CK-20-positive (1).

(104936-004)

M7019/EFG/KRM/2008.09.30 p. 1/4

Performance characteristics

Cells labelled by the antibody display a cytoplasmic staining pattern.

Normal tissues: The antibody labels normal urothelium (1, 4) and the mature epithelium lining the villi of duodenal mucosa (1). CK 20 immunoreactivity has not been detected in a number of non-epithelial tissues tested, such as smooth muscle, blood vessel walls, lymph nodes and tumour stroma (1).

Abnormal tissues: In colonic adenocarcinoma, the antibody labelled 26 of 27 tumours (96%). 21 (78%) displayed staining in more than 50% of the cells (3). Of 51 patients with urothelial papillomas of the bladder, 10 showed labelling of superficial urothelial cells (normal labelling) with the antibody. None of these patients showed further tumours during the follow up. By contrast, 30 patients (73%) showed CK 20 expression in all cell layers (abnormal labelling) and developed further tumours (4). 3 out of 4 metastatic gastric tumours were labelled by the antibody (5).

Focal immunostaining with the antibody (less than 10% positive cells) was observed in 9 of 65 lung, 1 of 20 breast, and 2 of 11 endometrial tumours. 1 of 19 clear-cell type renal cell carcinomas was positive in 10-50% of the cells (3). None of 10 primary ovarian carcinomas, and none of 5 mesotheliomas were found positive with the antibody (5).

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
References/ Références/ Literatur	<ol style="list-style-type: none">Moll R, Löwe A, Laufer J, Franke WW. Cytokeratin 20 in human carcinomas. A new histodiagnostic marker detected by monoclonal antibodies. Am J Pathol 1992;140:427-47. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. In: Herrmann, Harris, editors. Subcellular biochemistry. Volume 31, New York: Plenum Press; 1998. p. 205-62. Savera AT, Torres FX, Linden MD, Bacchi CE, Gown AM, Zarbo RJ. Primary versus metastatic pulmonary adenocarcinoma. An immunohistochemical study using villin and cytokeratins 7 and 20. Appl Immunohistochem 1996;4:86-94. Harnden P, Mahmood N, Southgate J. Expression of cytokeratin 20 redefines urothelial papillomas of the bladder. Lancet 1999; 353:974-77. Wauters CCAP, Smedts F, Gerrits LGM, Bosman FT, Ramaekers FCS. Keratins 7 and 20 as diagnostic markers of carcinomas metastatic to the ovary. Hum Pathol 1995;26:852-55. Moll R, Schiller DL, Franke WW. Identification of protein IT of the intestinal cytoskeleton as a novel type I cytkeratin with unusual properties and expression patterns. J Cell Biol 1990;111:567-80. Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. Cell 1982;31:11-24.
--	---

Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole				
REF	Catalogue number Référence du catalogue Bestellnummer	 <p>2 °C – 8 °C</p>	Temperature limitation Limites de température Zulässiger Temperaturbereich	 <p>Manufacturer Fabricant Hersteller</p>
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	LOT	Batch code Code du Lot Chargenbezeichnung	
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten		Use by Utiliser jusque Verwendbar bis	

(104936-004)	M7019/EFG/KRM/2008.09.30 p. 4/4
Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel. +45 44 85 95 00 Fax +45 44 85 95 95 CVR No. 33 21 13 17	