Monoclonal Mouse Anti-Human Follicular Dendritic Cell
 Clone CNA.42
Code No./ Code/ Code-Nr. M 7157
 Edition/ Edition/ Ausgabe 28.07.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Follicular Dendritic Cell, Clone CNA.42, is intended for use in immunocytochemistry. The antibody labels a 120 kDa antigen expressed on follicular dendritic cells. The antibody may be a useful tool for the identification of follicular dendritic cells, either forming network or scattered among malignant cells, in reactive lymphoid tissue, follicular lymphomas and Hodgkin's disease (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	Follicular dendritic reticulum cell (1).
Introduction	Follicular dendritic cells (FDCs) are present in primary lymphoid follicles and in the germinal centres of secondary lymphoid follicles of lymph nodes, spleen and tonsils. These cells are characterized by numerous cytoplasmic long branching processes extending between adjacent lymphoid cells that form a three-dimensional network (1). FDCs have important functions in the selection of memory B lymphocytes during germinal centre reactions. They present native antigens to potential memory cells, of which only B cells with high affinity B-cell receptors can bind. These B lymphocytes survive, whereas non-binding B cells undergo apoptotic cell death (2). The immunophysiological kinetics, developmental processes (including their cellular origin), morphological and functional transformation, terminal fate, and pathological behaviour of FDCs are still quite obscure (4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> CNA.42. <u>Isotype:</u> IgM, kappa. <u>Mouse IgG concentration:</u> 1 See label on vial.
Immunogen	Splenocytes from nude mice grafted with the human lymphoblastic CEM T-cell line (1).
Specificity	In Western blotting of lysed CEM cells under non-reducing conditions, the antibody labelled a band of 120 kDa. The CNA.42 antigen did not shift under reducing conditions. Neuraminidase digestion of tissue slides completely eliminated the reactivity of the antibody, suggesting a glycosylated epitope (1). As demonstrated by immunocytochemistry, the antibody cross-reacts with FDCs in bird, cat, cow, dog, goat, horse, monkey, mouse, rabbit, rat, sheep and swine (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308. Less optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, pH 9.0, code No. S 2368, and 10 mmol/L citrate buffer, pH 6.0, or pre-treatment of tissues with DakoCytomation Proteinase K, code No. S 3020. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Follicular Dendritic Cell, Clone CNA.42, may be used at a dilution range of 1:25 -1:50 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgM, code No. X 0942, diluted to the same mouse IgM concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. Follow the procedure enclosed with the selected visualization kit.
Performance characteristics	Cells labeled by the antibody display a staining pattern confined mainly to the cell membrane. <u>Normal tissues:</u> The antibody labelled FDCs in both primary and secondary lymphoid follicles in reactive lymph nodes, tonsils, white pulp of spleen and Peyer's patches. In thymus, the antibody labelled a subpopulation of large cortical thymocytes (1). In non-lymphoid organs the antibody labelled pancreatic islet cells, gastric chief cells, acini in salivary gland, Leydig cells of the testis, myelin sheaths, smooth muscles of the arterial wall, striated muscle fibers and endothelial cells (1). <u>Abnormal tissues:</u> In most of 289 cases of hematopoietic tumours of various types the antibody labelled either intact FDC networks or FDC networks dispersed among malignant cells (1).

The antibody labelled malignant cells in 6/129 (5%) cases of B-cell lymphomas, 30/184 (16%) cases of T-cell lymphomas, 23/105 (22%) cases of Hodgkin's disease, and 44/114 (39%) of diverse non-hematopoietic tumours such as different neurogenic tumours and carcinomas of gastrointestinal origin (1).

In 11/11 cases of Epstein-Barr virus related inflammatory pseudotumour-like FDC tumours, the antibody labelled FDCs, scattered mononuclear cells, mast cells as well as and some epithelial and mesenchymal cells (3).

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
1. Raymond I, Al Saati T, Tkaczuk J, Chittal S, Delsol G. CNA.42, a new monoclonal antibody directed against a fixative-resistant antigen of follicular dendritic reticulum cells. Am J Pathol 1997;151:1577-85.
2. van Nierop K, de Groot C. Human follicular dendritic cells: function, origin and development. Sem in Immunol 2002;14:251-7.
3. Cheuk W, Chan JKC, Shek TW, Chang JH, Tsou MH, Yuen NWF, et al. J. Inflammatory pseudotumor-like follicular dendritic cell tumor: a distinctive low-grade malignant intra-abdominal neoplasm with consistent Epstein-Barr virus association. Am J Surg Pathol 2001;25:721-31.
4. Maeda K, Matsuda M, Suzuki H, Saito H. Immunohistochemical recognition of human follicular dendritic cells (FDCs) in routinely processed paraffin sections. J Histochem Cytochem 2002;50:1475-85.

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