**Monoclonal Mouse
Anti-Human
Melanosome
Clone HMB45**

ENGLISH
Code M0634

Intended use

For In Vitro Diagnostic Use.

Monoclonal mouse anti-human melanosome, clone HMB45 (anti-melanosome, HMB45) is intended for laboratory use to identify qualitatively by light microscopy melanocytes with immature melanosome formation in normal and pathological paraffin-embedded tissue using immunohistochemical (IHC) test methods. Positive results aid in the classification of melanomas and melanocytic lesions and also aid in distinguishing metastatic amelanotic melanomas from other poorly differentiated tumors of uncertain origin. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Summary and explanation

The presence of the HMB45 antigen indicates active melanosome formation and thus melanocytic differentiation.⁵ It is also expressed in normal fetal melanocytes,^{2,3} but not in normal resting adult melanocytes, regardless of the degree of pigmentation.^{1,3,6} Upon activation, adult melanocytes can re-express the HMB45-defined antigen (as it is expressed in fetal melanocytes). Such melanocytes are activated by a variety of stimuli. For example, HMB45-positive cells have been detected in tissue overlying or adjacent to granulation tissue, hemangiomas, vessel-rich tumor stroma, and basal cell carcinoma.⁵⁻⁸ Hair follicles stain occasionally due to associated stimulated melanocytes.¹ Positive HMB45 staining has not been observed with melanocytes in lentigines or overlying fibroblastic proliferations such as keloids, dermatofibromas and old fibrotic hemangiomas.⁵ Non-melanocytic normal tissues do not react with the HMB45 antibody.

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: HMB45¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

M0634 may be used at a dilution of 1:50 when performing IHC using the EnVision™, EnVision™ Doublestain or LSAB™2 detection systems. These are guidelines only. Optimal antibody concentrations may vary depending on specimen and preparation method and should be determined by each individual laboratory.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Immunogen

Extract of pigmented melanoma metastases from lymph nodes

Specificity

Anti-melanosome, HMB45 has been shown to react with a 10 kD segment of a neuraminidase-sensitive sialylated glycoconjugate present in pre- and early-stage (immature) melanosomes.²⁻⁴

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions. Suggested diluent for IHC procedures:

Antibody Diluent (code S0809)

The following negative control is recommended for IHC procedures:

Mouse IgG₁ (code X0931)

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.⁹
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user.¹⁰ There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-melanosome, HMB45 can be used on formalin-fixed, paraffin-embedded tissue sections. Other fixatives that are compatible with anti-melanosome, HMB45 include methacarn, zinc formalin, neutral buffered formalin or B5. Pretreatment of tissue with proteolytic enzymes is not recommended.

Cryostat Sections and Cell Smears

Anti-melanosome, HMB45 can be used for labelling acetone-fixed cryostat sections or fixed cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-melanosome, HMB45 is cytoplasmic.

Product specific limitations

1. An HMB45-negative result does not exclude melanoma. Desmoplastic melanomas infrequently stain with the anti-melanosome, HMB45 antibody.^{4,11} Therefore it is recommended that this antibody be used in a panel of antibodies including anti-S-100. Anti-S-100 stains desmoplastic melanomas, but lacks the specificity of anti-melanosome, HMB45 for melanomas.
2. Anti-melanosome, HMB45 reacts with a melanocytic antigen, not a melanoma antigen. Anti-melanosome, HMB45 reactivity cannot differentiate between benign and malignant melanocytic lesions. Malignancy must be determined by other criteria.
3. Stimulated normal adult melanocytes (i.e. melanocytes overlying or adjacent to granulation tissue, hemangiomas, vessel-rich stroma and basal cell carcinomas) may express melanocytic antigens.⁵⁻⁸ In some basal cell carcinomas, rare HMB45-positive cells with dendritic morphology are scattered through-out and adjacent to the tumor. The origin of these HMB45-positive cells is unclear and may represent a "reactive" melanocytic population within the tumor.¹
4. Primary renal neoplasms are nonreactive, except for renal angiomyolipomas (RAML).^{12,13} An HMB45-related antigen located in the myeloid component of RAML appears to be similar to the melanocytic antigen in melanomas.¹²
5. Smooth muscle cells of pulmonary lymphangiomyomatosis (PLAM) exhibit an unusual phenotype characterized by anti-melanosome, HMB45 immunoreactivity distinct from other smooth muscle cell proliferations. Bonetti et al.¹⁴ tested 75 pathological lung transbronchial biopsies and only 6 PLAM specimens from 3 patients showed the presence of HMB45-positive cells.
6. Rare false-positive staining has been reported for nonmelanocytic tumors such as adenocarcinoma,^{15,16} olfactory neuroblastoma,¹² malignant lymphoma,¹⁶ clear cell tumor of the lung¹² and plasmacytoma¹⁷ with anti-melanosome, HMB45 ascites preparations.¹⁵⁻¹⁷ Nonspecific staining of normal breast epithelia, sweat glands, plasma cells and bronchial epithelia has also been observed with ascites-derived antibodies and appears as membranous and/or granular, cytoplasmic staining.¹⁷ These artifacts have been attributed to contaminated ascites preparation or inappropriately diluted antibody.^{4,15-17} No such staining has been reported with tissue culture supernatant preparations. Anti-melanosome, HMB45 is produced from tissue culture supernatant. If false-positive staining of nonmelanocytic neoplasia or normal epithelia occurs with anti-melanosome, HMB45, report the event to Dako's Technical Support Department.
7. This antibody does not work in Western blot assays.¹

Performance characteristics

Reproducibility

Anti-melanosome, HMB45 has been tested on serial sections of tissue specimens. Consistent staining results have been obtained with run-to-run and within-run antibody testing.

Normal Tissues

Normal adult tissues that exhibit positive staining with anti-melanosome, HMB45 include melanocytes (fetal and subset, melanocytes containing immature melanosomes), retinal pigment epithelia (prenatal and infantile). Negative tissues include adrenal gland, brain, breast, gallbladder, gastrointestinal tract, kidney, liver, lung, lymphoid tissue, mesenchymal cells, pancreas, peripheral nervous tissue, retinal pigment epithelia (adult), salivary gland, skin (melanocytes, normal resting, melanophages, Langerhans cells, keratinocytes, hair follicles, cutaneous nerves, sweat glands, sebaceous) and testis.

Abnormal Tissues

Anti-melanosome, HMB45 stains 245/256 (95.7%) of melanoma (excluding desmoplastic)^{1,3,5,6,11-14,18-20} and 245/291 (84.2%) of melanoma (including desmoplastic).^{1,3,5,6,11,18-20} Melanocytic atypical hyperplasia (2/2),¹ melanocytic neuroectoderm of infancy (1/1),⁶ renal angiomyolipoma (27/29)^{13,14} and various nevi (218/228)^{1,5,6,12} are stained by anti-melanoma, HMB45.

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Références

Literatur

- Gown AM, Vogel AM. Monoclonal antibodies to intermediate filament proteins of human cells: Unique and cross-reacting antibodies. J Cell Biol 1982;95(2 Pt 1):414-24
- Kapur RP, Bigler SA, Skelly M, Gown AM. Anti-melanoma monoclonal antibody HMB45 identifies an oncofetal glycoconjugate associated with immature melanosomes. J Histochem Cytochem 1992;40(2):207-12

3. Esclamado RM, Gown AM, Vogel AM. Unique proteins defined by monoclonal antibodies specific for human melanoma. *Amer J Surg* 1986;152(4):376-85
4. Taatjes DJ, Arendash-Durand B, von Turkovich M, Trainer TD. HMB-45 antibody demonstrates melanosome specificity by immunoelectron microscopy. *Arch Pathol Lab Med* 1993;117(3):264-8
5. Skelton HG, Smith KJ, Barrett TL, Lupton GP, Graham JH. HMB-45 staining in benign and malignant melanocytic lesions. *Amer J Dermatopathol* 1991;13(6):543-50
6. Colombari R, Bonetti F, Zamboni G, Scarpa A, Marino F, Tomezzoli A, Capelli P, Menestrina F, Chilosi M, Fiore-Donati L. Distribution of melanoma specific antibody (HMB-45) in benign and malignant melanocytic tumors. *Virch Arch A Pathol Anat* 1988;413(1):17-24
7. Smoller BR, Hsu A, Krueger J. HMB-45 monoclonal antibody recognizes an inducible and reversible melanocyte cytoplasmic protein. *J Cutan Pathol* 1991;18(5):315-22
8. Smoller BR, McNutt NS, Hsu A. HMB-45 recognizes stimulated melanocytes. *J Cutan Pathol* 1989;16(2):49-53
9. Department of Health, Education and Welfare, National Institute for Occupational Safety and Health, Rockville, MD. "Procedures for the decontamination of plumbing systems containing copper and/or lead azides." DHHS (NIOSH) Publ. No. 78-127, Current 13. August 16, 1976
10. Clinical Laboratory Improvement Amendments of 1988: Final Rule, 57FR7163. February 28, 1992
11. Wick MR, Swanson PE, Rocamora A. Recognition of malignant melanoma by monoclonal antibody HMB-45. An immunohistochemical study of 200 paraffin-embedded cutaneous tumors. *J Cutan Pathol* 1988;15(4):201-7
12. Ashfaq R, Weinberg AG, Albores-Saavedra J. Renal angiomyolipomas and HMB-45 reactivity. *Canc* 1993;71(10):3091-7
13. Pea M, Bonetti F, Zamboni G, Martignoni G, Riva M, Colombari R, Mombello A, Bonzanini M, Scarpa A, Ghimenton C, Donati LF. Melanocyte-marker-HMB-45 is regularly expressed in angiomyolipoma of the kidney. *Pathol* 1991;23(3):185-8
14. Bonetti F, Chiodera PL, Pea M, Martignoni G, Bosi F, Zamboni G, Mariuzzi GM. Transbronchial biopsy in lymphangiomyomatosis of the lung. HMB45 for diagnosis. *Amer J Surg Pathol* 1993;17(11):1092-102
15. Bonetti F, Pea M, Martignoni G, Mombello A, Colombari R, Zamboni G, Scarpa A, Piubello Q, Bacchi CE, Gown AM. False-positive immunostaining of normal epithelia and carcinomas with ascites fluid preparation of antimelanoma monoclonal antibody HMB45. *Amer J Clin Pathol* 1991;95(4):454-9
16. Bacchi CE, Gown AM. Letters to the Editor. Specificity of antibody HMB-45. *Arch Pathol Lab Med.* 1992;116(9):899-900
17. Leong AS-Y, Milios J. An assessment of a melanoma-specific antibody (HMB-45) and other immuno-histochemical markers of malignant melanoma in paraffin-embedded tissues. *Surg Pathol* 1989;2(2):137-45
18. Ordonez NG, Ji XL, Hickey RC. Comparison of HMB-45 monoclonal antibody and S-100 protein in the immunohistochemical diagnosis of melanoma. *Amer J Clin Pathol* 1988;90(4):385-90
19. Bacchi CE, Goldfogel GA, Greer BE, Gown AM. Paget's disease and melanoma of the vulva. *Gyn Oncol* 1992;46(2):216-21
20. Herman GE. Unpublished data. Dako 1992

REF Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	LOT Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595
www.dako.com

PT0039/Rev C

Edition 01/09