FLEX
Monoclonal Rabbit
Anti-Human AMACR
Clone 13H4
Ready-to-Use
(Link)

Code IR060

English

Intended use

For in vitro diagnostic use.

FLEX Monoclonal Rabbit Anti-Human AMACR, Clone 13H4, Ready-to-Use, (Link), is intended for use in immunohistochemistry together with Autostainer Link instruments. This antibody is useful for the identification of prostate adenocarcinomas (1-9). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Synonyms for antigen

Alpha methylacyl coenzyme A racemase, P504S (1, 2).

Summary and explanation

Monoclonal rabbit anti-AMACR (alpha methylacyl coenzyme A racemase) recognizes a 382 amino acid protein that was identified by cDNA library subtraction in conjunction with high throughput microarray screening of prostate adenocarcinoma (1). AMACR, also known as P504S, is an enzyme that is involved in bile acid biosynthesis and β -oxidation of branched-chain fatty acids (10). AMACR is expressed in cells of premalignant high grade prostatic intraepithelial neoplasia (HGPIN) and prostate adenocarcinoma, but is present at low or undetectable levels in glandular epithelial cells of normal prostate and benign prostatic hyperplasia (1-9). AMACR is expressed in normal non-prostatic epithelium and in carcinomas from outside the prostate (2, 11-23).

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Ready-to-use monoclonal rabbit antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L NaN_3 .

Clone: 13H4 (1-3)

Immunogen

Full length recombinant AMACR (1).

Specificity

The specificity of monoclonal rabbit anti-AMACR was evaluated by immunocytochemistry and Western blot analysis. Anti-AMACR positively bound formalin-fixed HEK293 cells which overexpressed AMACR, but was unreactive with cells transfected with an empty plasmid. In Western blots of lysates from primary prostate carcinoma samples, the antibody labels a 54 kDa protein corresponding to the expected molecular weight of AMACR (3).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN_3), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation including materials required but not supplied

The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm.

Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (10x) (Code K8000/K8004).

Paraffin-embedded sections: Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (10x) (Code K8000/K8004). Note: After staining, the sections must be dehydrated, cleared and mounted using permanent mounting medium.

Deparaffinized sections: Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.

Staining procedure including materials required but not supplied

The recommended visualization system is EnVision™ FLEX, High pH (Link) (Code K8000). The staining steps and incubation times are pre-programmed into the Autostainer Link software. The recommended reagent application volume is 1 x 200 µL or 2 x 150 µL per slide. Please refer to the proper Autostainer Link User Guide for detailed instructions on loading slides and reagents. If the protocols are not available on the used Autostainer platform, please contact Dako Technical Services. All incubation steps should be performed at room temperature.

Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, incubation times and the choice of EnVision™ FLEX/EnVision™ FLEX+ visualization system can be adjusted. A Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics."

Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin (Link) (Code K8008). Non-aqueous, permanent mounting medium is recommended.

Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include epithelial cells of the prostate glands and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Rabbit (Link) (Code IR600).

Staining interpretation

The cellular staining pattern is cytoplasmic.

Performance characteristics

Normal tissues:

Tissue Type (# Tested) (2-5, 11-14)	Positively Staining Tissue Elements
Adrenal (12) ²	0/12
Bladder (7) ²	0/7
Brain (5) ²	0/5
Breast (12) ²	0/12
Colon (24) ^{11,12}	20/24 Colonic surface epithelium, focal and luminal
Colon, hyperplastic polyp (28) ¹²	1/28
Endometrium (17) ^{2,11}	0/17
Gallbladder (10) ²	10/10 Epithelial cells
Heart (5) ²	0/5
Kidney (12) ²	12/12 Tubular epithelial cells
Liver (31) ^{2,13}	18/31 Hepatocytes
Lung (18) ^{2,11}	12/18 Bronchial epithelial cells
Lymph Node (7) ²	0/7
Ovary (14) ^{2,11}	0/14
Pancreas (15) ²	0/15
Prostate, normal (34) ^{2,11}	0/34
Prostate, benign (156) ^{3,4,5}	11/156
Salivary gland (7) ²	0/7
Skin (10) ^{2,11}	0/10
Small Intestine (18) ²	0/18
Spleen (12) ²	0/12
Stomach (13) ²	0/13
Stomach, non-neoplastic epithelium (44) ¹⁴	2/44
Testes (5) ²	0/5
Thyroid (5) ²	0/5

Abnormal tissues: Anti-AMACR is immunoreactive with the majority of prostate carcinomas tested (2-8,11).

Tissue Type (# Tested) (2, 11-23)	Positively Staining Tumors
Adrenal cortical tumor (20) ²	1/20
Basal cell carcinoma of skin (20) ²	0/20
Bile duct cholangiocarcinoma (14) ²	2/14
Breast adenocarcinoma (315) ^{2,11,15}	54/315
Carcinoid tumors; lung and GI (10) ²	0/10
Colon carcinomas (233) ^{11,12,16,17}	
Colon (35) ^{11,16}	16/35
Well differentiated (58) ¹²	45/58
Moderately differentiated (88) ¹²	66/88
Poorly differentiated (30) ¹²	11/30
Metastatic (22) ¹⁷	7/22
Endometrioid carcinoma (10) ²	0/10
Germ cell tumors (14) ²	0/14
Hepatocellular carcinoma (72) ^{2,13}	59/72
Lung carcinoma (300) ^{2,11,18}	
Carcinoma (28) ¹¹	4/28
Adenocarcinoma (151) ^{2,18}	75/151
Squamous cell carcinoma (121) ¹⁸	27/121
Melanoma (41) ^{2,11}	2/41
Neuroendocrine carcinoma; GI, lung and liver (167) ^{2,18}	95/167
Ovary adenocarcinoma (76) ^{2,11,17}	2/76
Pancreas adenocarcinoma (13) ²	1/13
Pleural mesothelioma (16) ²	0/16
Prostate adenocarcinoma (640) ^{2-8,11}	587/640
Renal tumors (105) ^{19,20}	
Clear cell carcinoma (77) ^{19,20}	17/77
Chromophobe (24) ^{19,20}	1/24
Oncocytoma (29) ^{19,20}	4/29
Bellini/urothelial (5) ¹⁹	0/5
Mucinous and spindle cell carcinoma (5) ¹⁹	5/5
Papillary (96) ^{19,20}	96/96
Sarcomatoid (15) ²⁰	0/15
Salivary gland tumors (28) ²	1/28
Small cell carcinoma; lung and skin (15) ²	0/15
Soft tissue epithelioid sarcoma (12) ²	1/12
Soft tissue synovial sarcoma (6) ²	0/6
Squamous cell carcinoma; skin and mucosa (25) ²	0/25
Stomach adenocarcinoma (413) ^{2,14,21,22}	262/413
Undifferentiated carcinoma (27) ²	0/27
Thymoma (8) ²	0/8
Thyroid tumors (54) ²	0/54
Urothelial carcinoma (50) ^{2,20,23}	
Urothelial carcinoma (17) ²⁰	2/17
Transitional cell carcinoma, bladder (29) ²	9/29
Clear cell adenocarcinoma (4) ²³	4/4

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

**References
Bibliographie
Literaturhinweise**

1. Xu J, Stolk JA, Zhang X, Silva SJ, Houghton RL, Matsumura M, Vedvick TS, Leslie KB, Badaro R, Reed SG. Identification of differentially expressed genes in human prostate cancer using subtraction and microarray. *Canc Res*. 2000; 60:1677
2. Jiang Z, Fanger GR, Woda BA, Banner BF, Algate P, Dresser K, Xu J, Chu PG. Expression of α -methylacyl-CoA racemase (P504S) in various malignant neoplasms and normal tissues: A study of 761 cases. *Human Pathology*. 2003; 34(8):792
3. Jiang Z, Woda BA, Rock KL, Xu Y, Savas L, Khan A, Pihan G, Cai F, Babcook JS, Rathanaswami P, Reed SG, Xu J, Fanger GR. P504S: a new molecular marker for the detection of prostate carcinoma. *Am J Surg Pathol*. 2001; 25(11):1397
4. Jiang Z, Wu C, Woda BA, Dresser K, Xu J, Fanger GR, Yang XJ. P504S/ α -methylacyl-CoA racemase: a useful marker for diagnosis of small foci of prostatic carcinoma on needle biopsy. *Am J Surg Pathol*. 2002; 26(9):1169
5. Yang XJ, Wu C, Woda BA, Dresser K, Tretiakova M, Fanger GR, Jiang Z. Expression of α -methylacyl-CoA racemase (P504S) in atypical adenomatous hyperplasia of the prostate. *Am J Surg Pathol*. 2002; 26(7):921
6. Zhou M, Jiang Z, Epstein JI. Expression and diagnostic utility of alpha-methylacyl-CoA-racemase (P504S) in foamy gland and pseudohyperplastic prostate cancer. *Am J Surg Pathol*. 2003; 27(6):772
7. Beach R, Gown AM, de Peralta-Venturina MN, Folpe AL, Yaziji H, Salles PG, Grignon DJ, Fanger GR, Amin MG. P504S immunohistochemical detection in 405 prostatic specimens including 376 18-gauge needle biopsies. *Am J Surg Pathol*. 2002; 26(12):158
8. Wu C-L, Yang XJ, Tretiakova M, Patton KT, Halpern EF, Woda BA, Young RH, Jiang Z. Analysis of α -methylacyl-coA racemase (P504S) expression in high-grade prostatic intraepithelial neoplasia. *Hum Pathol* 2004;35(8):1008-13
9. Jiang Z, Woda B, Wu C-L, Yang X. Discovery and clinical application of a novel prostate cancer marker. *Am J Clin Pathol* 2004;122:275-89
10. Ferdinandusse S, Denis S, Ijst L, Dacremont G, Waterham HR, Wanders JA. Subcellular localization and physiological role of α -methylacyl-CoA racemase. *Journal of Lipid Research*. 2000; 41:1890
11. Nassar A, Amin MB, Sexton DG, Cohen C. Utility of α -methylacyl coenzyme A racemase (P504S) antibody) as a diagnostic immunohistochemical marker for cancer. *Appl Immunohistochem Mol Morphol* 2005;13(3):252-5
12. Jiang Z, Fanger GR, Banner BF, Woda BA, Algate P, Dresser K, Xu J, Reed SG, Rock KL, Chu PG. A dietary enzyme: alpha-methylacyl-CoA racemase/P504S is overexpressed in colon carcinoma. *Cancer Detect Prev*. 2003;27(6):422
13. Li W, Cagle PT, Botero RC, Liang JJ, Zhang Z, Tan D. Significance of overexpression of alpha methylacyl-coenzyme A racemase in hepatocellular carcinoma. *J Exp Clin Cancer Res* 2008;27:21
14. Lee WA. α -methylacyl-CoA-racemase expression in adenocarcinoma, dysplasia and non-neoplastic epithelium of the stomach. *Oncology* 2006;71:246-50
15. Witkiewicz AK, Varambally S, Shen R, Mehra R, Sabel MS, Ghosh D, Chinnaiyan AM, Rubin MA, Kleer CG. α -methylacyl CoA racemase protein expression is associated with the degree of differentiation in breast cancer using quantitative image analysis. *Canc Epidemiol Biomarkers Prev* 2005;14(6):1418-23
16. Osunkoya AO, Netto GJ, Epstein JI. Colorectal adenocarcinoma involving the prostate: report of 9 cases. *Human Pathol* 2007;38:1836-41
17. Logani S, Oliva E, Arnell PM, Amin MB, Young RH. Use of novel immunohistochemical markers expressed in colonic adenocarcinoma to distinguish primary ovarian tumors from metastatic colorectal carcinoma. *Mod Pathol* 2005;18:19-25
18. Shilo K, Dracheva T, Mani H, Fukuoka J, Sesterhenn IA, Chu W-S, Shih JH, Travis WD, Franks TJ. α -methylacyl CoA racemase in pulmonary adenocarcinoma, squamous cell carcinoma, and neuroendocrine tumors. *Arch Pathol Lab Med* 2007;131:1555-60
19. Molinié V, Balaton A, Rotman S, Mansouri D, de Pinieux I, Homsy T, Guillou L. Alpha-methyl CoA racemase expression in renal cell carcinomas. *Human Pathol* 2006;37:698-703
20. Tretiakova MS, Sahoo S, Takahashi M, Turkyilmaz M, Vogelzang NJ, Lin F, Krausz T, Teh BT, Yang XJ. Expression of alpha-methylacyl-coA racemase in papillary renal cell carcinoma. *Am J Surg Pathol* 2004;28(1):69-76
21. Truong CD, Li W, Feng W, Cagle P, Khoury T, Sadir A, Xie K, Yao J, Tan D. Alpha-methylacyl-CoA racemase expression is upregulated in gastric adenocarcinoma: a study of 249 cases. *Int J Clin Exp Path* 2008;1:518-23
22. Cho EY, Kim K-M, Park CK, Kim JJ, Sohn TS, Kim DW. AMACR is highly expressed in gastric adenomas and intestinal-type carcinomas. *APMIS* 2007;115:713-8
23. Sun K, Huan Y, Unger PD. Clear cell adenocarcinoma of urinary bladder and urethra: another urinary tract lesion immunoreactive for P504S. *Arch Pathol Lab Med* 2008;132:1417-22

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Contains sufficient for <N> tests Contenu suffisant pour <n> tests Inhalt ausreichend für „n“ Ansätze	Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA
 Use by Utiliser jusque Verwendbar bis	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	 EC REP Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel 805 566 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763 Tel +45 4485 9500 Fax +45 4485 9595 www.dako.com

PT0020/ Rev C