Monoclonal Mouse Anti-Human CD138, Clone MI15

Code No./ Code/ Code Nr. **C 7256** **APC-Conjugated**
 Code No./ Code/ Code Nr. **R 7229** **RPE-Conjugated**

ENGLISH

Intended use For in vitro diagnostic use.
 C 7256 and R 7229 are intended for use in flow cytometry. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

Synonym for antigen Syndecan-1 (1).

Introduction CD138 is a transmembrane proteoglycan, bearing heparan sulphates and sometimes chondroitin sulfates. It consists of a 30.5 kDa core protein and five glycosaminoglycans. CD138 is a member of the syndecan family. The reported Mr of CD138 ranges from 85 000 to 92 000, presumably depending on the glycosaminoglycan pattern of the molecule (1).
 Within the haematopoietic system, CD138 is mainly confined to late stages of B-cell differentiation. It is typically expressed at high cellular density in normal and malignant plasma cells (2) and in lymphoplasmacytoid cells (3). In mature tissue CD138 is expressed by simple and stratified epithelia, fibroblasts, stratified keratinocytes and endothelial cells (1).
 The function of CD138 is largely unknown. However, the intracellular tail of CD138 seems to bind together the cytoskeleton with cellular components (1), and the extracellular part of the molecule seems to bind different ligands. Presumably, CD138 plays a role in cellular functions, such as proliferation, programmed cell death, cell-matrix and cell-cell adhesion (2), for example it mediates adhesion of myeloma cells to type I collagen (4).
 CD138 expression is reduced during malignant transformation of various epithelia (5), and the antigen is rapidly lost (shed) by myeloma cells entering into apoptosis, thus CD138 is a marker of viable myeloma cells (1, 4). CD138 is present in malignancies involving terminally differentiated plasma cells, and it is associated with multiple myeloma (4, 6, 7). Together with BCL-6 and IRF4, CD138 may be of value for the diagnosis of AIDS-related non-Hodgkin lymphomas and HIV-related Hodgkin lymphomas (8). Surface-expressed CD138 has been reported to be sensitive to proteolytic cleavage (1).

Reagent provided The Anti-CD138 conjugates, C 7256 and R 7229, have been produced from a purified monoclonal mouse antibody. The conjugates are provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN₃, pH 7.2. Each vial contains 100 tests (10 µL of conjugate for up to 5 x 10⁵ U266 cells).
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code No.	Fluorochrome	Negative Control Code No.
R 7229	RPE (R-Phycoerythrin)	X 0928
C 7256	APC (Allophycocyanin)	X 0968

Immunogen A mixture of U266 and XG-1 human myeloma cell lines (6, 9).

Specificity Anti-CD138, MI15, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Kobe 1996), and studies by a number of different laboratories confirmed its reactivity with CD138 (1).
 The epitope recognized by anti-CD138, MI15, is found within the ectodomain (extracellular part) of the CD138 core protein. Blocking experiments have revealed that the epitope recognized by MI15 seems to be partly overlapping the epitope recognized by another antibody to CD138, clone B-B4 (1, 2). Anti-CD138, MI15, reacts strongly with CD138 expressed by multiple myeloma-derived cell lines, such as U266 and XG-1 (9).
In blood: No significant reaction is seen with leukocytes from normal peripheral human blood (< 5%).
In bone marrow: No significant reaction is seen in normal bone marrow (< 5%). CD138 is absent from CD34 positive stem cells (1). In bone marrow from multiple myeloma patients, all plasma cell types are labelled, including reticular, polymorphous, asynchronous and blastic plasma cells (6).

Precautions 1. For professional users.
 2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
 3. This product contains material of animal origin, therefore the product should be handled as potentially infectious.

Storage Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

Staining procedure

1. Collect venous blood into a test tube containing an anticoagulant.
2. Isolate mononuclear cells by centrifugation on a separation medium. Alternatively, lyse the red cells after step 6.
3. Wash the mononuclear cells twice with RPMI 1640 or PBS, pH 7.2-7.4.
4. Mix 100 µL cell suspension with 10 µL fluorochrome-conjugated Anti-CD138.
5. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
6. Incubate in the dark at 4 °C for 30 minutes.
7. Wash twice with PBS containing 2% BSA. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.

8. Analyse on a flow cytometer.

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

Product-specific limitations

Owing to the sensitivity of CD138 towards proteolytic cleavage, the processing of samples, such as freezing or storage before staining with the antibody, may play an important role in the outcome of the staining (1).

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
1. Wijdenes J, Clément C, Klein B, Dore J-M. New B-cell CD antigens. BC29: CD138 (syndecan-1) workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 249-52.
2. Gattei V, Godeas C, Degan M, Rossi FM, Aldinucci D, Pinto A. Characterization of anti-CD138 monoclonal antibodies as tools for investigating the molecular polymorphism of syndecan-1 in human lymphoma cells. Br J Haematol 1999;104:152-62.
3. Sebestyén A, Berczi L, Mihalik R, Paku S, Matolcsy A, Kopper L. Syndecan-1 (CD138) expression in human non-Hodgkin lymphomas. Br J Haematol 1999;104:412-9.
4. Jourdan M, Ferlin M, Legouffe E, Horvathova M, Liautard J, Rossi JF, et al. The myeloma cell antigen syndecan-1 is lost by apoptotic myeloma cells. Br J Haematol 1998;100:637-46.
5. Inki P, Jalkanen M. The role of syndecan-1 in malignancies [review]. Ann Med 1996;28:63-7.
6. Costes V, Magen V, Legouffe E, Durand L, Baldet P, Rossi J-F, et al. The Mi15 monoclonal antibody (anti-syndecan-1) is a reliable marker for quantifying plasma cells in paraffin-embedded bone marrow biopsy specimens. Hum Pathol 1999;30:1405-11.
7. van Zaanen HCT, Vet RJWM, de Jong CM, von dem Borne AEG, van Oers MHJ. A simple and sensitive method for determining plasma cell isotype and monoclonality in bone marrow using flow cytometry. Br J Haematol 1995;91:55-9.
8. Carbone A, Gloghini A, Larocca LM, Capello F, Pierconti F, Canzonieri V, et al. Expression profile of MUM1/IRF4, BCL-6, and CD138/syndecan-1 defines novel histogenetic subsets of human immunodeficiency virus-related lymphomas. Blood 2001; 97:744-51.
9. Horvathova M, Gaillard JP, Liautard J, Duperray C, Lavabre-Bertrand T, Bourquard P, et al. B30.6. Identification of novel and specific antigens of human plasma cells by mAb. In: Schlossman SF, Boumsell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995, volume two. p. 713-4.

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