

Anti-CCNB3 Product Datasheet
Antibody

PRODUCT SPECIFICATION

Certificate of Analysis

Product Name	Anti-CCNB3
Product Number	HPA000496
Lot Number	A82690
Gene Description	cyclin B3
Clonality	Polyclonal
Isotype	IgG
Host	Rabbit
Antigen Sequence	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: KMCASQRKQSCQEESLAVQDVNMEEDSFFMESMSFKKKPKTEESIPTHKL SSLKKKCTIYGKICHFRKPPVLQTTICGAMSSIKKPTTEKETLFFQELSVL QEKHTTEHEMS
Purification Method	Affinity purified using the PrEST antigen as affinity ligand
Verified Species	Human
Reactivity	
Recommended Applications	IHC (Immunohistochemistry) - Antibody dilution: 1:200 - 1:500 - Retrieval method: HIER pH6 ICC-IF (Immunofluorescence) - Fixation/Permeabilization: PFA/Triton X-100 - Working concentration: 1-4 µg/ml
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Unit size	100 µl
Concentration	0.3 mg/ml
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.
Shipping	Normally shipped at ambient temperature
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user. See Protocols under Technical Aids www.sigma.com/tech aids .

Certificate is system generated and does not require a signature

 Product of Sweden. MSDS available
For research use only. Not for drug, household or other uses