

Anti-CA6 Product Datasheet  
Antibody

## Certificate of Analysis

### PRODUCT SPECIFICATION

Product Name	Anti-CA6
Product Number	HPA028692
Lot Number	A105057
Gene Description	carbonic anhydrase VI
Clonality	Polyclonal
Isotype	IgG
Host	Rabbit
Antigen Sequence	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: SEISGSEHTVDGIRHVIEIHIVHYNKYKSYDIAQDAPDGLAVLAAFVEV KNYPENTYYSNFISHLANIKYPGQRTTLTGLDVQDMLPRNLQHYYTY
Purification Method	Affinity purified using the PrEST antigen as affinity ligand
Verified Species	Human
Reactivity	
Recommended Applications	IHC (Immunohistochemistry) - Antibody dilution: 1:20 - 1:50 - Retrieval method: HIER pH6
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Unit size	100 µl
Concentration	0.05 mg/ml
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.
Shipping	Normally shipped at ambient temperature
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user. See Protocols under Technical Aids <a href="http://www.sigma.com/techaids">www.sigma.com/techaids</a> .

Certificate is system generated and does not require a signature

 Product of Sweden. MSDS available  
For research use only. Not for drug, household or other uses

**SIGMA-ALDRICH™**

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich Publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.