

2019-03-29

9th Fl., No.108, Jhouzih St., Neihu District, Taipei City, Taiwan, R.O.C.

Tel: + 886-2-8751-1888 Fax: + 886-2-6602-1218

E-Mail: sales@abnova.com

377.1-PRAGO00001_FedEx_CZ_

Certificate of Analysis

Product Name : MTAP monoclonal antibody (M01), clone 2G4**Catalog Number :** H00004507-M01**Size :** 0.1MG**Lot Number :** J2211-2G4**Immunogen :** MTAP (AAH18625, 1 a.a. ~ 154 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.**Sequence :** MASGTTTTAVKIGIIGGTGLDDPEILEGRTEKYVDTPFGKPSDALILGKIKNVDCILLARHGRO
HTIMPSKVNYQANIWALKEEGCTHVIVTTACGSLREEIQPGDIVIIDQFIDSHVRRACFPFTFH
HDCFQRPPQKPSRSHCVSCATCRTMS**Clonality :** Monoclonal**Host :** Mouse**Reactivity :** Human,Mouse,Rat**Concentration :** 0.5mg/ml**Entrez Gene ID :** 4507**Gene Symbol :** MTAP**Gene Alias :** MSAP|c86fus**Isotype :** IgG1 Kappa**Storage Buffer :** In 1x PBS, pH 7.4**Storage Instruction :** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.**Applications :** ELISA,IF,S-ELISA,WB-Ce,WB-Re
(See application key <http://www.abnova.com/products/seekey.htm>)**Images :** See www.abnova.com for image details by key in **Catalog Number** under **Abnova Search**.
(Image may not be available for certain products)Approved by : 

Quality Assurance Representative

If you find any defect of this product on receipt, please inform us within 10 days, we will replace the faulty product.

Certificate is system generated and does not require a signature.

Product of Taiwan, MSDS available.

For R&D use only Not for drug, household or other uses. The Product is INEDIBLE and CAN NOT be used as food additive.

Abnova (Taiwan) Corporation warrants that its products conform to the information contained in this and other Abnova (Taiwan) Corporation Publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.