

Safety Data Sheet (SDS) According to the REACH Regulation (EC) No. 1907/2006

Issuing Date: 2017-07-10

Revision Date: 2023-10-03

Version: 2

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product No 38313
Product name GPNMB (E4D7P) XP® Rabbit mAb

Contains

Chemical name	Index No.	CAS No
glycerol (30-60)	Not Listed	56-81-5
sodium azide (<0.02)	011-004-00-7	26628-22-8

1.2. Relevant identified uses of the substance or mixture and uses advised against

Identified uses For Research Use Only. Not for Use in Diagnostic Procedures.

1.3. Details of the supplier of the safety data sheet

Importer	Manufacturer
Cell Signaling Technology Europe B.V. Dellaertweg 9b 2316 WZ Leiden The Netherlands TEL: +31 (0)71 7200 200 FAX: +31 (0)71 891 0019	Cell Signaling Technology, Inc. 3 Trask Lane Danvers, MA 01923 United States TEL: +1 978 867 2300 FAX: +1 978 867 2400

Website www.cellsignal.com
E-mail Address info@cellsignal.eu

1.4. Emergency telephone number

CHEMTREC 24 hours a day, 7 days a week, 365 days a year
+1 703 527 3887 (INTERNATIONAL) +1 800 424 9300 (NORTH AMERICA)

Europe 112

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Regulation (EC) No. 1272/2008

This substance is classified as not hazardous according to regulation (EC) 1272/2008 [CLP]

2.2. Label elements

Signal word

None.

Hazard statement(s)

None.

Precautionary statement(s)

None.

2.3. Other hazards

May produce an allergic reaction.

*For the full text of the H-phrases & EUH-phrases mentioned in this Section, see Section 16***SECTION 3: Composition/information on ingredients****Chemical nature** Mixture

Chemical name	CAS No	Weight-%	EC No	Classification (1272/2008)	REACH Registration Number
glycerol	56-81-5	30-60	200-289-5	-	no data available
sodium azide	26628-22-8	<0.02	247-852-1	Acute Tox. 2 (H300) Aquatic Acute 1 (H400) Aquatic Chronic 1 (H410) (EUH032)	no data available

*For the full text of the H-phrases & EUH-phrases mentioned in this Section, see Section 16***SECTION 4: First aid measures****4.1. Description of first aid measures****General advice**

Use first aid treatment according to the nature of the injury. When symptoms persist or in all cases of doubt seek medical advice.

Inhalation

Move to fresh air.

Skin contact

Wash skin with soap and water.

Eye contact

Rinse thoroughly with plenty of water, also under the eyelids.

Ingestion

Clean mouth with water and afterwards drink plenty of water.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain, or flushing.

4.3. Indication of any immediate medical attention and special treatment needed**Notes to physician**

Treat symptomatically.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable Extinguishing Media**

Use extinguishing measures that are appropriate to local circumstances and the

Unsuitable Extinguishing Media surrounding environment.
No information available.

5.2. Special hazards arising from the substance or mixture

No information available.

5.3. Advice for firefighters

Wear self-contained breathing apparatus and protective suit. Use personal protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

For non-emergency personnel Avoid contact with skin, eyes and clothing. Use personal protective equipment. For personal protection see section 8.
For emergency responders Use personal protection recommended in Section 8.

6.2. Environmental precautions

Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. Prevent entry into waterways, sewers, basements or confined areas.

6.3. Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.
Methods for cleaning up Soak up with inert absorbent material. Pick up and transfer to properly labeled containers.

6.4. Reference to other sections

See Sections 8 & 13 for additional information.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Handle in accordance with good industrial hygiene and safety practice. Wear personal protective equipment. Avoid contact with skin, eyes and clothing. Remove and wash contaminated clothing before re-use.

7.2. Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

7.3. Specific end use(s)

Use as a laboratory reagent.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Chemical name	European Union	United Kingdom	France	Spain	Germany
glycerol		STEL 30 mg/m ³ TWA 10 mg/m ³	TWA 10 mg/m ³	TWA 10 mg/m ³	Ceiling / Peak: 400 mg/m ³ TWA: 200 mg/m ³

sodium azide	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ S*	STEL 0.3 mg/m ³ TWA 0.1 mg/m ³ Skin	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ P*	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ S*	TWA: 0.2 mg/m ³ Ceiling / Peak: 0.4 mg/m ³
Chemical name	Italy	Portugal	Netherlands	Finland	Denmark
glycerol		TWA 10 mg/m ³		TWA 20 mg/m ³	
sodium azide	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ Pelle*	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ Ceiling 0.29 mg/m ³ Ceiling 0.11 ppm C(A4) P*	Huid* STEL 0.3 mg/m ³ TWA 0.1 mg/m ³	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ iho*	TWA 0.1 mg/m ³ H*
Chemical name	Austria	Switzerland	Poland	Norway	Ireland
glycerol		SS-C** TWA 50 mg/m ³ STEL 100 mg/m ³	TWA 10 mg/m ³		TWA 10 mg/m ³ STEL 30 mg/m ³
sodium azide	H* STEL 0.3 mg/m ³ TWA 0.1 mg/m ³	TWA 0.2 mg/m ³ STEL 0.4 mg/m ³	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³	TWA 0.1 mg/m ³ STEL 0.1 mg/m ³	TWA 0.1 mg/m ³ STEL 0.3 mg/m ³ Skin

8.2. Exposure controls

Appropriate engineering controls

Showers, eyewash stations, and ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection

Safety glasses with side-shields

Skin protection

Wear protective gloves and protective clothing

Hand protection

Impervious gloves.

Other

Wear suitable protective clothing.

Respiratory protection

When workers are facing concentrations above the exposure limit they must use appropriate certified respirators.

Environmental Exposure Controls

No information available.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	Liquid - Clear
Color	Colorless
Odor	No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	7.5	@ 20 °C
Melting point/freezing point	No information available	No information available
Boiling point or initial boiling point and boiling range	No information available	No information available
Flash point	No information available	No information available.
Evaporation rate	No information available	No information available
Flammability	No information available	No information available
Upper/lower flammability or explosive limits	No information available	No information available
Vapor pressure	No information available	No information available
Relative vapor density	No information available	No information available
Density and/or relative density	No information available	No information available
Solubility	No information available.	No information available
Partition coefficient: n-octanol/water	No information available	No information available

Autoignition temperature	No information available	No information available
Decomposition temperature	No information available	No information available.
Viscosity	No information available	No information available
Explosive properties	No information available	No information available
Oxidizing properties	No information available	No information available

9.2. Other information

Softening point	No information available
Molecular Weight	No information available
Solubility in other solvents	No information available
VOC content	No information available
Liquid Density	No information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No information available.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization	Hazardous polymerization does not occur.
Hazardous reactions	None under normal processing.

10.4. Conditions to avoid

Extremes of temperature and direct sunlight. Over a period of time, sodium azide may react with copper, lead, brass, or solder in plumbing systems to form an accumulation of the HIGHLY EXPLOSIVE compounds of lead azide & copper azide.

10.5. Incompatible materials

Strong oxidizing agents, Strong acids.

10.6. Hazardous decomposition products

Nitrogen oxides (NO_x).

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

This product is for experimental uses only. The product has not been completely analyzed and all of the hazards may not be known. Please use caution while handling this product.

Chemical name	LD50 Oral	LD50 Dermal	LC50 Inhalation
glycerol	= 12600 mg/kg (Rat)	> 10 g/kg (Rabbit)	> 570 mg/m ³ (Rat) 1 h
sodium azide	= 27 mg/kg (Rat)	= 20 mg/kg (Rabbit) = 50 mg/kg (Rat)	-

Information on likely routes of exposure

Inhalation Avoid breathing vapors or mists.
Eye contact Avoid contact with eyes.
Skin contact Avoid contact with skin.
Ingestion Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.

Symptoms Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain, or flushing.

Skin corrosion/irritation No information available.
Serious eye damage/eye irritation No information available.
Sensitization No information available.
Mutagenic effects No information available.
Carcinogenicity No information available.

Reproductive toxicity No information available.
STOT - single exposure No information available.
STOT - repeated exposure No information available.
Aspiration Hazard No information available.

11.2. Information on other hazards

No information available.

SECTION 12: Ecological information

12.1. Toxicity

Chemical name	Toxicity to algae	Toxicity to fish	Toxicity to daphnia and other aquatic invertebrates
glycerol	-	LC50 51 - 57 mL/L (Oncorhynchus mykiss) 96 h	EC50 500 mg/L (Daphnia magna) 24 h
sodium azide	EC50 0.35 mg/L (Pseudokirchneriella subcapitata) 96 h	LC50 0.8 mg/L (Oncorhynchus mykiss) 96 h LC50 5.46 mg/L (Pimephales promelas) 96 h LC50 0.7 mg/L (Lepomis macrochirus) 96 h	LC100 1 mg/L (Orconectes rusticus) 96 h

12.2. Persistence and degradability

No information available

12.3. Bioaccumulative potential

Bioaccumulation

Chemical name	Octanol-Water Partition Coefficient
glycerol	-1.76

Bioconcentration factor (BCF) No information available.

12.4. Mobility in soil

No information available.

12.5. Results of PBT and vPvB assessment

No information available.

12.6. Endocrine disrupting properties

This product does not contain any known or suspected endocrine disruptors

12.7. Other adverse effects

No information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste from residues / unused products	Dispose of in accordance with local regulations.
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal.
Other information	Waste codes should be assigned by the user based on the application for which the product was used.

SECTION 14: Transport information

IMDG/IMO

14.1 UN number	Not regulated
14.2 UN proper shipping name	Not regulated
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	None
14.6 Special precautions for user	None
14.7 Maritime transport in bulk according to IMO instruments	Not regulated

ADR/RID

14.1 UN number	Not regulated
14.2 UN proper shipping name	Not regulated
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	None
14.6 Special precautions for user	None

IATA

14.1 UN number	Not regulated
14.2 UN proper shipping name	Not regulated
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	None
14.6 Special precautions for user	None

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Candidate List of Substances of Very High Concern for Authorization Information

This product does not contain Substances of Very High Concern (SVHC).

SEVESO Directive Information

This product does not contain substances identified in the SEVESO Directive.

International inventories

TSCA 8(b)	Complies
DSL/NDSL	Complies
EINECS/ELINCS	Complies
ENCS	-
IECSC	Complies
KECL	-
PICCS	-
AICS	Complies

International inventories legend

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory
DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List
EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances
ENCS - Japan Existing and New Chemical Substances
IECSC - China Inventory of Existing Chemical Substances
KECL - Korean Existing and Evaluated Chemical Substances
PICCS - Philippines Inventory of Chemicals and Chemical Substances
AICS - Australian Inventory of Chemical Substances

15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out

SECTION 16: Other information

Full text of H-Statements referred to under Sections 2 and 3

H300 - Fatal if swallowed
H400 - Very toxic to aquatic life
H410 - Very toxic to aquatic life with long lasting effects
EUH032 - Contact with acids liberates very toxic gas

Classification procedure: Expert judgment and weight of evidence determination.
Issuing Date: 2017-07-10
Revision Date: 2023-10-03

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.