

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 3.0 Revision Date 28.10.2009

Print Date 29.10.2010

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

Product name : Anti-LCT, antibody produced in rabbit

Product Number : HPA007408
Brand : Sigma

Company : Sigma-Aldrich spol. s r.o.
Sokolovska 100/94
CZ-186 00 PRAHA 8

Telephone : +420246003200
Fax : +420246003292
Emergency Phone # : +420 603 441 855
E-mail address : eurtechserv@sial.com

2. HAZARDS IDENTIFICATION

Not a hazardous substance or mixture according to EC-directives 67/548/EEC or 1999/45/EC.

3. COMPOSITION/INFORMATION ON INGREDIENTS

CAS-No.	EC-No.	Index-No.	Classification	Concentration
Glycerol				
56-81-5	200-289-5	-	-	40 %
Sodium dihydrogenorthophosphate				
7558-80-7	231-449-2	-	-	1,08 %
Sodium chloride				
7647-14-5	231-598-3	-	-	0,51 %
Water				
7732-18-5	231-791-2	-	-	58,391 %
Sodium azide				
26628-22-8	247-852-1	011-004-00-7	T+, N, R28 - R32 - R50/53	0,02 %

For the full text of the R-phrases mentioned in this Section, see Section 16.

4. FIRST AID MEASURES

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

5. FIRE-FIGHTING MEASURES**Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self contained breathing apparatus for fire fighting if necessary.

6. ACCIDENTAL RELEASE MEASURES**Personal precautions**

Avoid breathing vapors, mist or gas.

Environmental precautions

Do not let product enter drains.

Methods for cleaning up

Keep in suitable, closed containers for disposal.

7. HANDLING AND STORAGE**Handling**

Normal measures for preventive fire protection.

Storage

Store in cool place. Keep container tightly closed in a dry and well-ventilated place.

Recommended storage temperature: -20 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION**Personal protective equipment****Respiratory protection**

Respiratory protection is not required. Where protection is desired, use multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection

For prolonged or repeated contact use protective gloves.

Eye protection

Safety glasses

Hygiene measures

General industrial hygiene practice.

9. PHYSICAL AND CHEMICAL PROPERTIES**Appearance**

Form liquid

Safety data

pH no data available

Melting point no data available

Boiling point no data available

Flash point	no data available
Ignition temperature	no data available
Lower explosion limit	no data available
Upper explosion limit	no data available
Water solubility	no data available

10. STABILITY AND REACTIVITY

Storage stability

Stable under recommended storage conditions.

Materials to avoid

Strong bases, Strong oxidizing agents

11. TOXICOLOGICAL INFORMATION

Acute toxicity

no data available

Irritation and corrosion

no data available

Sensitisation

no data available

Chronic exposure

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Signs and Symptoms of Exposure

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

Potential Health Effects

Inhalation	May be harmful if inhaled. May cause respiratory tract irritation.
Skin	May be harmful if absorbed through skin. May cause skin irritation.
Eyes	May cause eye irritation.
Ingestion	May be harmful if swallowed.
Target Organs	Kidney,

12. ECOLOGICAL INFORMATION

Elimination information (persistence and degradability)

no data available

Ecotoxicity effects

no data available

Further information on ecology

no data available

13. DISPOSAL CONSIDERATIONS

Product

Observe all federal, state, and local environmental regulations.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION**ADR/RID**

Not dangerous goods

IMDG

Not dangerous goods

IATA

Not dangerous goods

15. REGULATORY INFORMATION**Labelling according to EC Directives**

Further information:

The product does not need to be labelled in accordance with EC directives or respective national laws.

16. OTHER INFORMATION**Text of R-phrases mentioned in Section 3**

R28	Very toxic if swallowed.
R32	Contact with acids liberates very toxic gas.
R50/53	Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Further information

Copyright 2009 Sigma-Aldrich Co. License granted to make unlimited paper copies for internal use only. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Co., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.