

BG-8 Monoclonal Antibody (SIGNET)

Catalog Number: SIG-3317

Notes:

Positive tissue (human): Colon adenocarcinoma

Storage:

Store between 2-8°C.

References:

Yaziji H, Battifora H, Barry TS, Hwang HC, Bacchi CE, McIntosh MW, Kussick SJ, Gown AM. Evaluation of 12 antibodies for distinguishing epithelioid mesothelioma from adenocarcinoma: identification of a three-antibody immunohistochemical panel with maximal sensitivity and specificity. *Mod Pathol* 19(4):514-23, 2006.

Ordóñez NG. The diagnostic utility of immunohistochemistry and electron microscopy in distinguishing between peritoneal mesotheliomas and serous carcinomas: a comparative study. *Mod Pathol* 19(1):34-48, 2006.

Ordóñez NG. The diagnostic utility of immunohistochemistry in distinguishing between epithelioid mesotheliomas and squamous carcinomas of the lung: a comparative study. *Mod Pathol* 19(3):417-28, 2006.

Ordóñez NG. Value of thyroid transcription factor-1, E-cadherin, BG8, WT1, and CD44S immunostaining in distinguishing epithelial pleural mesothelioma from pulmonary and nonpulmonary adenocarcinoma. *Am J Surg Pathol* 24(4): 598-606, 2000.

Zhu ZM, Kojima N, Stroud MR, et al. Monoclonal antibody directed to Le(y) oligosaccharide inhibits implantation in the mouse. *Biol Reprod* 52:903-12, 1995.

Gaffney RA, Schaeffer AJ, Duncan JL. Lewis blood group antigen expression by cultured normal ureteral epithelial cells. *J Urol* 148:1341-1346, 1992.

Jordon D, Jagirdar J, Kaneko M. Blood group antigens, Lewisx and Lewisy in the diagnostic discrimination of malignant mesothelioma versus adenocarcinoma. *1494-1494 Pathol* 135(5):931-7, 1989.

Warranty/Conditions:

Covance products may not be resold or modified for resale without prior written approval. All sales are subject to Covance Antibody Products Terms and Conditions of Sale.

Related Products:

- Ultra Streptavidin (USA) Horseradish peroxidase (HRP) 50 Test, Multi-Species, DAB Detection Kit (SIGNET), Catalog Number: SIG-32250
- BG-5 Monoclonal Antibody (SIGNET), Catalog Number: SIG-3314
- BG-6 Monoclonal Antibody (SIGNET), Catalog Number: SIG-3315

Antibody Products: +1.800.922.2226
+1.510.594.3800
Email: ab.products@covance.com

The logo for Covance, featuring the word "COVANCE" in a bold, blue, sans-serif font. A red swoosh underline is positioned beneath the letters "VANCE".