




SAB4502272 SIGMA

# Anti-BCOR antibody produced in rabbit

~1 mg/mL, affinity isolated antibody, buffered aqueous solution  
 Synonym: BCL-6 corepressor, BCoR, KIAA1575

SDS



## Properties

Related Categories	<a href="#">Alphabetical Index</a> , <a href="#">Antibodies</a> , <a href="#">B-BI</a> , <a href="#">Primary Antibodies</a>
species reactivity	mouse, human
application(s)	immunohistochemistry: suitable western blot: suitable
clone	polyclonal
concentration	~1 mg/mL
antibody form	affinity isolated antibody
form	buffered aqueous solution
mol wt	antigen mol wt 192 kDa
shipped in	wet ice
storage temp.	-20°C
Gene Information	human ... <a href="#">BCOR(54880)</a>
biological source	rabbit
conjugate	unconjugated
NCBI accession no.	<a href="#">Q6W2J9</a>

## Description

### Immunogen

Synthesized peptide derived from human BCOR

### General description

Anti-BCOR Antibody detects endogenous levels of total BCOR protein.

### Physical form

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Price and Availability

SKU-Pack Size	Availability	Price (EUR/CZK)	Quantity
SAB4502272-100UG	Only 1 left in stock (more on the way) - FROM	9,760.50	<input type="text" value="0"/>

[Bulk orders?](#)

