

Anti-AFF2 Product Datasheet
Antibody

PRODUCT SPECIFICATION

Certificate of Analysis

Product Name	Anti-AFF2
Product Number	HPA003139
Lot Number	A08752
Gene Description	AF4/FMR2 family, member 2
Clonality	Polyclonal
Isotype	IgG
Host	Rabbit
Antigen Sequence	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: SSTSDSNTDQEETLQIKVLPPCIISGNTAKSKEICGASLTLSTLMSSG SNNLSISNEEPTFSPIPVMQTEILSPLRDHENLKNLWVKIDL DLLSRVP GHSSLHAAPAKPDHKETATKPKRQTAVTAVEKPAPKGKRKH
Purification Method	Affinity purified using the PrEST antigen as affinity ligand
Verified Species Reactivity	Human
Recommended Applications	IHC (Immunohistochemistry) - Antibody dilution: 1:20 - 1:50 - Retrieval method: HIER pH6
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Concentration	0.05 mg/ml
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user. See Protocols under Technical Aids www.sigma.com/techaids .

Certificate is system generated and does not require a signature

 Product of Sweden. MSDS available
For research use only. Not for drug, household or other uses

SIGMA-ALDRICH™

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich Publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.