




BOB.1 (SP92)



Catalog Number:	760-4593
Ordering Code:	06433308001
Quantity:	50 tests
Controls:	Tonsil
Clone Name:	SP92
Species:	Rabbit Monoclonal
Localization:	Nuclear, Cytoplasmic
Regulatory Status:	IVD

This antibody is intended for in vitro diagnostic (IVD) use. BOB.1 (SP92) antibody is intended for qualified laboratories to qualitatively identify by light microscopy the presence of associated antigens in sections of formalin-fixed, paraffin-embedded tissue sections using IHC test methods. Use of this antibody is indicated as an aid in the identification of B-cells in germinal centers, mantle cells, and differentiation of lymphomas within the context of an antibody panel, the patient's clinical history, and other diagnostic tests evaluated by a qualified pathologist.

© 2016 Ventana Medical Systems, Inc. [Terms and Conditions](#) [Privacy Policy](#) 15.04.2014

This may be a regulated product; it may not be available in all geographies. Please consult your local representative for availability and restrictions. This website contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any legal process, regulation, registration or usage in the country of your origin.